

Dodge

LOOKING FOR A **FORMALDEHYDE-FREE CHOICE** FOR YOUR EMBALMING? **DODGE HAS IT.**

Freedom Art: Preserves, disinfects, deodorizes • Pleasantly scented • Non-spotting, uniform diffusion
Imparts vibrant, natural tissue color • Superior clearing of discolorations • Non-firming

Freedom Cav: No gas buildup • No liquid buildup • No odors
Pleasant apple scent • Preservation that lasts



For Professional
Embalming Use Only

For optimum results, the following guidelines are recommended minimums.

200 lbs. or less

1 container of **Freedom Art**
3 bottles of **ProfFlow**
1 bottle of **Rectifiant**
Enough water to make a 2
gallon solution.

200 lbs. and heavier

Begin with the standard suggested mixture.
For each additional 50 lbs. of weight, add:
24 oz. **Freedom Art** plus
one 16 oz. bottle of **Rectifiant**, and
one 16 oz. bottle of **ProfFlow** to enough
water to make a 1 gallon solution.

Follow the arterial embalming with cavity treatment using a minimum of 2 bottles
of **Freedom Cav**. Distribute 1 bottle into the thoracic region and 1 bottle into the
abdominal region.

Freedom Art
Case of 4 (120 oz. Jug)
Catalog No. 108011

120 Oz. Jug
Catalog No. 108012

ProfFlow
Case of 24
Catalog No. 340040

Rectifiant
Case of 24
Catalog No. 320010

Freedom Cav
Case of 24
Catalog No. 241776